


# Two Arm Wrist Rest

## SECTION VI PARTS LIST



MA319701

| Item | Part No.       | Description                                               | Qty. | Item | Part No.       | Description                                                | Qty. |
|------|----------------|-----------------------------------------------------------|------|------|----------------|------------------------------------------------------------|------|
|      | 9A222001       | Two Arm Wrist Rest Assembly<br>(Includes Items 1 thru 37) | 1    | 20   | • 040-0250-88  | • Screw (Apply Loctite #042-0025-00)                       | 4    |
| 1    | • 530-1059-00  | • Wrist Tube Weldment                                     | 1    | 21   | • 029-1100-04  | • One Arm Headlock Assembly<br>(Includes Items 22 thru 36) | 1    |
| 2    | • 553-0093-02  | • Round Tube Insert                                       | 2    | 22   | • 040-0250-11  | • Screw                                                    | 1    |
| 3    | • 040-0250-120 | • Screw (Apply Loctite #042-0025-00)                      | 2    | 23   | • 040-0010-14  | • Set Screw                                                | 2    |
| 4    | • 057-0355-00  | • Pivot Spacer                                            | 2    | 24   | • 030-0131-07  | • Headlock Head                                            | 1    |
| 5    | • 040-0250-119 | • Set Screw                                               | 2    | 25   | • 042-0002-08  | • Roll Pin                                                 | 1    |
| 6    | • 040-0250-118 | • Screw                                                   | 2    | 26   | • 057-0048-01  | • Pivot Bar                                                | 1    |
| 7    | • 551-0706-00  | • Pivot Stop Block                                        | 1    | 27   | • 051-0080-100 | • Draw Braw                                                | 1    |
| 8    | • 529-0689-00  | • Wrist Rest Post Assembly                                | 1    | 28   | • 030-0123-03  | • Head Pivot Handle Weldment                               | 1    |
| 9    | • 529-0688-00  | • Knob Assembly                                           | 1    | 29   | • 042-0002-04  | • Roll Pin                                                 | 2    |
| 10   | • 530-1056-00  | • Clamp Weldment                                          | 1    | 30   | • 041-0250-12  | • Jam Nut                                                  | 2    |
| 11   | • 520-0055-06  | • Handle (Apply Lube #064-0001-00)                        | 1    | 31   | • 040-0250-70  | • Set Screw                                                | 2    |
| 12   | • 030-1032-00  | • Head Plate Weldment                                     | 1    | 32   | • 040-0250-11  | • Screw                                                    | 1    |
| 13   | • 028-0431-00  | • Headrest Cushion (Specify Color)                        | 1    | 33   | • 040-0250-139 | • Screw (Apply Loctite #042-0025-00)                       | 4    |
| 14   | • 040-0010-102 | • Screw (Apply Loctite #042-0025-00)                      | 4    | 34   | • 030-0132-03  | • Headlock Weldment                                        | 1    |
| 15   | • 040-0010-92  | • Screw (Apply Loctite #042-0025-00)                      | 3    | 35   | • 057-0048-00  | • Pivot Bar                                                | 1    |
| 16   | • 030-0124-03  | • Headboard Pivot Assembly                                | 1    | 36   | • 051-0080-00  | • Draw Braw                                                | 1    |
| 17   | • 041-0010-07  | • Acorn Nut                                               | 4    | 37   | • 045-0000-01  | • Stainless Washer                                         | 2    |
| 18   | • 029-0103-02  | • Headboard Pivot Assembly                                | 1    | 38   | • 003-0796-00  | • Installation Instructions                                | 1    |
| 19   | • 042-0001-00  | • Roll Pin                                                | 1    |      |                |                                                            |      |

Always Specify Model &amp; Serial Number